
Autopsy Table ST 10/350 Fixed Mounted – BASIC *line*





The table is set on an installation column and two height-adjustable feet to to equalize out any unevenness. The autopsy table is completely made of stainless steel (material no. 1.4301- ANSI 304). The work surface features a high, profiled surround and a seamlessly welded basin, $400 \times 500 \times 200$, with a tube overflow are installed to the end of the table. The table top has a negative inclination arranged diagonally to the drain. The drains from the supporting surface and the basin are brought together in one drain connection. All power connections are built into the base and can be easily accessed by means of an inspection opening.

The Autopsy table ST 10/350 model includes the following features:

- · Knee-operated mixing tap for cold and warm water
- 3 m shower hose with hand sprinkler
- 2 splash proof electrical outlets 230 Volt German style
- The safety requirements for drinking water as per DIN 1988 part 4 are to be provided by the customer
- · Electrical and water connections inside the installation column are to be provided by the customer

Technical Data:

Dimensions W/D/H: $2600 \times 850 \times 850 \text{ mm}$ Voltage: 230 Volt / 1 Ph / 50 Hz Carring capacity: max. 200 kg

Type: ST 10/350
Article No.: 010.010.350

KUGEL medical GmbH & Co. KG Hermann-Köhl-Str. 2A

DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0 Fax +49/941/20 86 48-29 www.KUGEL-medical.de





Autopsy Table ST 10/350 Fixed Mounted - BASIC line







Ph: 0800 456 217 info@tekmed.co.nz

www.tekmed.co.nz

TekMed NZ TekEquipment

Ph: 1300 368 138 www.tekequipment.com.au

info@tekequipment.com.au



