Stainless Steel Body Tray LM 80/100



The body tray serves to accommodate, to store and to transport bodies of deceased. A seamless deepdrawn, extremely rugged structure of ground stainless steel (material no. 1.4301 -ANSI 304) without welded corners and edges on the inside ensures perfect cleaning. Two closed skids at the bottom of the body tray serves for noiseless transport on the roller tracks of the lift and transport truck or on the shelf system, and prevents moreover lateral tipping. A big sliding handle is mounted at each front end of the body tray. The body tray is absolutely resistant to cleaning agents and disinfectants.

Technical Data:

Dimensions W/D/H: Carrying capacity: Dead weight: Material: Type: Article No.: 2080 x 620 x 90 mm max. 200 kg approx. 22 kg Stainless steel (Material No. 1.4301 – ANSI 304) LM 80/100 080.010.010

Additional models see back page

KUGEL medical GmbH & Co. KG Hermann-Köhl-Str. 2A DE – 93049 Regensburg / Germany Phone +49/941/208648-0 Fax +49/941/208648-29 www.KUGEL-medical.de





Subject to technical modifications as required by products improvements (Version 01/2014)

Additional Models

Body tray LM 80/100 - 750-250kg

Dimensions W/D/H: Carrying capacity: Dead weight: Material: Typ: Article No.: 2080 x 750 x 90 mm max. 250 kg approx. 30 kg Stainless steel (Material No 1.4301 – ANSI 304) LM 80/100-750-250kg 080.010.009

Body tray LM 80/100 - 750-350kg

Dimensions W/D/H: Carrying capacity: Dead weight: Material: Typ: Article No.: 2080 x 750 x 90 mm max. 350 kg approx. 30 kg Stainless steel (Material No 1.4301 – ANSI 304) LM 80/100-750-350kg 080.010.082

Body tray LM 80/100-800-350 kg

Dimensions W/D/H: Carrying capacity: Dead weight: Material: Type: Article No.: 2080 x 800 x 90 mm max. 350 kg approx. 50 kg Stainless steel (Material No. 1.4301 – ANSI 304) LM 80/100-800-350kg 089.100.002

Body tray LM 80/101 - sideway loading

Dimensions W/D/H: Carrying capacity: Dead weight: Material: Type: Article No.: 1900 x 620 x 100 mm max. 200 kg approx. 22 kg Stainless steel (Material No. 1.4301 – ANSI 304) LM 80/101 080.010.014

Option:

Drain for Body Tray, Article No. 080.010.030

TekEquipment

Ph: 1300 368 138

www.tekequipment.com.au

info@tekequipment.com.au

TekMed NZ[®]

Ph: 0800 456 217 info@tekmed.co.nz www.tekmed.co.nz

KUGEL medical GmbH & Co. KG Hermann-Köhl-Str. 2A DE – 93049 Regensburg / Germany Phone +49/941/208648-0 Fax +49/941/208648-29 www.KUGEL-medical.de





Subject to technical modifications as required by products improvements (Version 01/2014)