

Staining Table FT 850-AK-ULS - Standing Workplace



This specialized work table provides a pollutant-free laboratory environment while staining histological slide preparations. Airflow travels over the entire working area with an integrated fan as well as a permanganate carbon filter and a HEPA filter. The air is vacuumed out downwards which provides optimal removal of unsanitary vapours. MAC (maximum allowable concentration) values are below target throughout the entire vacuum-effected area. Stable stainless steel construction with perforated working plate (with and without openings for glass cuvettes). Sub construction with height-adjustable feet. 3-side Plexiglas cover with stainless steel backside. Soundproofed silent exhaust fan with operation hour counter. Removable water container approx. 0.8 Liter.

Technical Data:

Dimensions W/D/H :	850 x 850 x 1750 mm
Air speed:	min. 0,4 m/sec on the perforated working plate
Material:	Stainless steel (material no. 1.4301 – ANSI 304)
Waste drain:	DN 50 (optional)
Voltage:	230 Volt / 1 Ph / 50
Type:	FT-850-AK-ULS
Article No.:	060.260.150

More information see back page

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■


Staining Table FT 850-AK-ULS - Standing Workplace



Perforated working area with openings for glass cuvettes (optional without openings)

Accessories:

Water tap with long hose
Type: EHMB-300
Article no. 010.010.040



Carbon filter (spare part)

Type: AKF-500-500/200
Article No. 060.260.026

TekEquipment™

Ph: 1300 368 138

www.tekequipment.com.au

info@tekequipment.com.au

TekMed NZ™

Ph: 0800 456 217

info@tekmed.co.nz

www.tekmed.co.nz

HEPA Filter (spare part)

Type: HEPA H13-305-610
Article no. 060.265.020

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical